



ELECTRODES & ACCESSORIES

Maximum result both in terms of stability and precision, also under the most difficult conditions. The line of electrodes is extremely complete and makes it possible to meet any monitoring requirement.

- For adults or for paediatric use
- With solid gel or liquid gel
- For normal monitoring or for stress tests
- FOAM or TNT/CLOTH or MICROPORE backing
- Round- shaped with different diameters, oval, rectangular or square. This line of accessories offers complementary or similar articles that are in great demand. Also in the case of accessories, the continuous research and exchange of information with the final user have allowed us to manufacture high quality and extremely reliable products.

TOP TRACE LINE

MONITORING & RESTING ECG'S



NM 50 RFI

Ø 50 mm.	Dimensions	Ø 55 mm.
Foam	Backing	Foam
Solid	Gel	Liquid
Ag/AgCl	Sensor	Ag/AgCl



NM 55 RLI

40 pcs / pouch - 800 pcs / box - 9600 pcs / carton

Packaging

25 pcs / pouch - 500 pcs / box - 6000 pcs / carton

The most typical sizes. The most demanded supports. A proved quality always useful.



NM 42 DFI

45x42 mm.	Dimensions	42x36 mm.
Foam	Backing	Foam
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl



NM 3642 OFI

50 pcs / pouch - 1200 pcs / box - 14400 pcs / carton

Packaging

50 pcs / pouch - 1200 pcs / box - 14400 pcs / carton

Two multifunctional models, suitable for adult and pediatric ECG's.



NP 30 RFI

Ø 30 mm.	Dimensions	28x44 mm.
Foam	Backing	Foam
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl



NM 2844 RFS

30 pcs / pouch - 900 pcs / box - 10800 pcs / carton

Packaging

50 pcs / pouch - 1250 pcs / box - 15000 pcs / master carton

The classical and largely used electrode for pediatric ECG's.
Excellent quality ensured also on "critical" patients.



NM 3351 OFI

51x33 mm.	Dimensions	28x44 mm.
Foam	Backing	Foam
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl - Banana plug



NM 2844 FBI

30 pcs / pouch - 900 pcs / box - 10800 pcs / carton

Packaging

30 pcs / pouch - 600 pcs / box - 7200 pcs / carton

TOP TRACE LINE

STRESS TEST



ST 50 RFI

Ø 50 mm.	Dimensions	Ø 50 mm.
Foam	Backing	Foam
Solid	Gel	Liquid
Ag/AgCl	Sensor	Ag/AgCl

40pcs/pouch-800pcs/box- 9600pcs/carton

Packaging

ST 50 RLI



Ø 50 mm.

Foam	Backing	Foam
Solid	Gel	Liquid
Ag/AgCl	Sensor	Ag/AgCl

30 pcs / pouch - 600 pcs / box - 7200 pcs / carton

The most suitable electrodes for ECG's under difficult conditions.
Easy to apply and use • Strong adhesiveness • Do not irritate the skin • Highly conductive gel

MRI NON - MAGNETIC



RT 50 RFI

Ø 50 mm.	Dimensions	23x34 mm.
Foam	Backing	Silicone coated pet film
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl

30pcs/pouch-600pcs/box- 7200pcs/carton

Packaging

ROUTINE ECG MONITORING

MEDTAB



23x34 mm.

Foam	Backing	Silicone coated pet film
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl

100 pcs / pouch - 500 pcs / box - 5000 pcs / carton

This electrode entirely replaces suction chests and clamps to patient's advantage both in terms of hygiene and health. To be attached to alligator clip connectors.

MONITORING / PRE-WIRED ELECTRODES



RE 30 LWF

30 mm.
female 4 mm

Dimensions

Foam	Backing	Foam
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl

3 pcs / pouch - 300 pcs / box

Available also with 22 mm diameter Model RE 22 LWF (with 4 mm female socket)

Packaging

RE 22 LWS

22 mm.
safety 1,55 mm

Dimensions

Foam	Backing	Foam
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl

3 pcs / pouch - 300 pcs / box

Available also with 30 mm diameter Model RE 30 LWS (with 1,55 mm safety plug)

Packaging



RT 30 LWF

30 mm.
female 4 mm

Dimensions

Foam	Backing	Foam
Solid	Gel	Solid
Carbon	Sensor	Carbon

3 pcs / pouch - 300 pcs / box

Available also with 22 mm diameter Model RT 22 LWF (with 4 mm female socket)

Packaging



RT 22 LWS

22 mm.
safety 1,55 mm

Dimensions

Foam	Backing	Foam
Solid	Gel	Solid
Carbon	Sensor	Carbon

3 pcs / pouch - 300 pcs / box

Available also with 30 mm diameter Model RT 30 LWS (with 1,55 mm safety plug)

Packaging

MEDI SENSOR LINE

MONITORING & RESTING ECG'S



MODEL F 50-S

Ø 50 mm.	Dimensions	Ø 26 mm.
Foam	Backing	Foam
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl



F 26-S

30 pcs / pouch – 1200 pcs / carton **Packaging** 50 pcs / pouch – 2000 pcs / carton
 Available also with liquid gel Model F 50-L

HOLTER / LONG MONITORING



MODEL W 601-S

Ø 55 mm.	Dimensions	Ø 55 mm.
TNT	Backing	MICROPORE
Solid	Gel	Solid
Ag/AgCl	Sensor	Ag/AgCl



MODEL T 601-S

30 pcs / pouch – 1500 pcs / carton **Packaging** 30 pcs / pouch – 1500 pcs / carton
 Available also with Ø 50 mm Model W 601 and also with liquid gel Model W 601-L

STRESS TEST



MODEL FS 50-L

Ø 50 mm.	Dimensions	Ø 55 mm.
Foam	Backing	Foam
Liquid	Gel	Liquid
Ag/AgCl	Sensor	Ag/AgCl



MODEL F 55-L

30 pcs / pouch – 1200 pcs / carton **Packaging** 30 pcs / pouch – 1200 pcs / carton

MULTIFUNCTION

MODEL F 50-S BLUE



MODEL FS 501-S

Ø 35x45 mm	Dimensions	Ø 50 mm
Foam	Backing	Foam
Solid	Gel	Solid
Ag/AgCl	Sensor	Carbon



50 pcs / pouch – 2000 pcs / carton **Packaging** 30 pcs / pouch – 1200 pcs / carton

Ideal for all ECG monitoring, radiolucent application included!
 AVAILABLE ALSO FOR PEDIATRIC, Ø 26 mm, MODEL F 26-S BLUE

SUCTION CHEST ELECTRODES



Ag/AgCl sensor to guarantee maximum reliability in signal transmission.

- Available in three sizes: Ø15 mm for pediatric use / Ø 24 mm and Ø 30 mm for adult use.

RUBBER LIMB & THORACIC STRAPS



For attaching electrodes • for limb use - one hole (15 mm x 40 cm / 25 mm x 45 cm), - two holes (32mm x 45cm) - for thoracic use - one hole (32mm x 135cm), six holes (100mm x 135cm).

CLAMP ELECTRODES



Classical-designed clamps with an improved quality. Ag/AgCl to ensure excellent signal transmission • Silver component to prevent distortions and to increase efficiency • Universal connection.

ELECTRO-SURGICAL PLATES



ELECTROSURGICAL PLATES WITH A SUPERIOR QUALITY

• Made of FOAM suitable to obtain perfect results with all instruments • manufactured and tested in compliance with EN 60601-1/EN 60601-2/AAMI-HF 18 • available with one zone and two zones and REM alarm for adults (adhesive surface 215 sq cm), for children (surface 129 sq cm) and for new born (surface 63 sq cm). **All the models are available even equipped with cables.**

DEFIBRILLATOR PADS



A wide range of pads suitable for all the defibrillators available. Raw material quality and attention to details resulted into an absolutely reliable range of products. Universal without cable or equipped with cable, having a specific connector to comply with all instruments.



ELECTRODES FOR TENS AND EMG



• Superior quality • available with single or double cable • lead wire with 2 mm socket or press-stud connection • manufactured in compliance with EN ISO 9001 and AAMI standards.

PULSE OXIMETER SENSORS



A broad range of sensors to monitor and detect the arterial oxygen saturation. Excellent quality and proven accuracy are the characteristics that better describe a whole series of disposable and reusable sensors. Suitable for all the main fittings spread on the market.

ELECTRODES & ACCESSORIES FOR EEG & EMG



- A) BRIDGE ELECTRODES** Made of silver/silver chloride to optimize signal transmission. Slightly curved support to improve head adaption.
- B) SILICONE HEADCAPS** With 2 or 4 adjustable bands to fit for various requirements. Anatomically-designed earmuffs. Available in three sizes (Small-Medium-Large);
- C) SPECIAL INK FOR EEG** Suitable & excellent for all devices (Siemens excluded), allowing a perfect trace without blurs. Black colour- 1 litre bottle
- D) DISC ELECTRODES** Available in silver/silver chloride or gold - 10 mm diameter- With 120-150-200-250 cm lead wire & DIN 42802 safety plug - 10 electrodes per set with coloured lead wires (10 colours available);
- E) DISPOSABLE CONCENTRIC NEEDLES FOR ELECTROMYOGRAPHY** Ergonomic handgrip and extremely sharp point - Equipped with wings for automatic unfastening -Shielded, light & flexible lead wire with 5 pole DIN 45322 connector and entirely shielded handgrip. Different diameters (from 0.30 to 0.60 mm) and different lengths (from 25 to 75 mm) - 25 pcs per packing.

MOUTHPIECES



Made of material particularly resistant to saliva or manipulation, perfectly hygienic, without residues and non-allergic • can be used on most spirometers of the most important brands present on the market • available in the external diameters 22-26-28-30 and 34 mm, but in case of large orders the customer can request a specific diameter • boxes of 500 pcs • an extremely handy dispenser is available on customer's demand.

